



Technical data



Equipment and accessories

standard ● optional ○

	MTS2-E	MTS2-H	MTS3-E	MTS3-H	TTS3-E	TTS3-H	TTS4	TTS8	BTS	BTS2
durable, soft foam upholstery (8 cm) – SER-SFU8-S	●	○	●	○	●	○	●	●	●	●
durable, firm foam upholstery (8 cm) – FFU8-S	○	○	○	○	○	○	○	○	○	○
breathing hole with welded and taped seams (plug included)	●	○	●	○	●	○	○	○	○	○
synchronously adjustable three-part head section – TVL-S	○	○	○	○	○	○	○	●	/	/
paper roll holder under head section (not available with TVL-S) – NZP1-S	○	○	○	○	○	○	○	/	/	/
paper roll holder at table base – NZP2-S	○	○	○	○	○	○	○	○	/	/
PVC cover for leg section – PVC8-S	○	○	○	○	○	○	○	○	○	○
retractable castors with soft release (Ø 75 mm) – SER-RC5	●	○	●	○	●	○	●	○	●	○
hand switch	/	/	/	/	/	/	●	●	/	/
circular switch	○	/	○	/	○	/	○	○	/	/
safety switch – EVV-S	○	/	○	/	○	/	○	○	○	○
battery backup – BB-S	○	/	○	/	○	/	○	○	○	○
plug connector for potential equalization – PPEQ	○	○	○	○	○	○	○	○	○	○
safe working load 250 kg	●	○	●	○	●	○	●	○	●	○



in f Bring treatment to a higher level!

Novak M d.o.o.

Pod smrekami 5  
SI - 1218 Komenda

Email: info@novak-m.si  
Tel: +386 1 5620117  
Fax: +386 1 5620118

www.novak-m.com



massage and physiotherapy  
treatments



reddot winner 2020



# MASSAGE AND THERAPY TABLES S

**SIMPLICITY** on one hand and **RICHNESS** on the other are the two distinguishable features of the renewed S-line design.

**Upgraded functionality** and continuation of the awarded design make massage and therapy tables the **perfect support** for different **massages**, **physiotherapy treatments** and **therapies**. The new features simplify the use and make tables even more **comfortable** for both patients and practitioners.



## So the cleaning becomes easier

To ensure **simple cleaning** and **high hygiene standards**, the **UPHOLSTERY** is **SEAMLESS** and made of **high-quality BIOCOMPATIBLE imitation leather**. **SEAMS** around the breathing hole are **WELDED** and **TAPED** to **prevent penetration** of any **liquids, massage lotions** and **oils** in particular. The **ABS plastic COVER** that protects the table base is **screwless** and thus completely **smooth**. The head and leg section **ADJUSTMENT LEVERS** with a mechanism **integrated** in the frame are **easier to press** and **clean**.

## To maximize the comfort and ease the use

To maximize the patient comfort, the **8 cm thick PADDING** is made of **two-layer foam**. It is **durable and comfortable** along the **whole lying surface** which is now even wider due to the **padding finished at a right angle**. To fully support the treatments, you can choose between **soft** and **firm** foam padding. The **BREATHING HOLE** with **memory foam** on the edges ensures perfect comfort in prone position. A **matching PLUG** is standard too. Optionally, tables can include a **THREE-PART HEAD SECTION** that allows **practitioners** an **unhindered access** when they treat patients seated at the head section. With the **NEW MECHANISM** for armrest adjustment, only a **push by one hand** is enough to lower or raise both armrests at once. **Lowering** and **raising** massage and therapy tables by the new **CIRCULAR SWITCH** is as simple as it can be. From any side of the table, one can adjust height by **pressing** the switch **down** or **pushing** it **up** completely effortless. If one chooses, the height on **Massage tables S2, S3** and on **Therapy table S3** can be adjusted hydraulically by pressing the **FOOT PEDAL** on either side of the table. The standard **SAFE WORKING LOAD** has been increased to **250 kg**, which contributes to the **perfect stability**. Practitioners can rely on the table even if they apply strong pressure during treatments. Due to the **NEW SOFTWARE** that enables **synchronous movement** of the middle two lying sections, they can **position** their **patients easier** and **quicker**. The **HAND SWITCH** includes two buttons for the **prestored** completely **flat** and **Trendelenburg positions** and a **memory function** to store a position of practitioner's choice. The table stands firmly on its **RUBBER FEET**. To move the table across the room, it features **RETRACTABLE CASTORS**. Due to the **shock absorber**, the release of the castors is softened.

## For various treatments, massages and therapies

The **BREATHING HOLE** with a **MATCHING PLUG** and **ELECTRIC HEIGHT ADJUSTMENT** are now the **standard equipment** of the renewed S-line massage and therapy tables. Two massage table models offer a **wide and comfortable lying surface** to perform massages and therapies. If treatments demand positions with negative section incline to allow the **Jackknife** or **Trendelenburg position**, 3- , 4- and 8-section therapy table models are available. The **fully equipped** Therapy table S8 standardly includes a **three-part head section** and **sturdier foldable side wings** integrated in the back section. Bobath table offers **WIDE SURFACE** for rehabilitation treatments. Now both, one and two-section model, include a **HAND SWITCH** to adjust the height and with a two-section Bobath the head section as well.